

TransferMan®

Product features of TransferMan

- Routine manipulator for the micromanipulation of suspension cells
- Easy to operate
- Dynamic control of manipulator movement, with two selectable speeds
- Two independent storage positions
- Home function ensures easy capillary exchange




- Easy-to-use keypad

Product features of TwinTip-Holder

- For special manipulation techniques (e.g. biopsy)
- Fits on all TransferMan and TransferMan NK workstations
- For the independent control of two capillaries on one TransferMan
- Special mechanics (patent pending) for extremely rapid and safe exchange between the two capillaries



- Workstation with Leica Microscope



Ordering information

Article	Order no.
TransferMan NK Proportional micromanipulator for suspension cells and microdissection	5178 000.014
TransferMan Dynamic micromanipulator for suspension cells	5177 000.010
CellTram Oil Manual piston pump for precise control of cell transfer (sperm, ES cells, etc.). Consists of CellTram Oil, 1 m pressure tube, capillary holder for capillaries with 1 mm outer diameter, 25 ml paraffin oil.	5176 000.025
CellTram Air Manual piston pump for safe holding of suspended cells (e.g. oocytes). Consists of CellTram Air, 1 m pressure tube, capillary holder for capillaries with 1 mm outer diameter.	5176 000.017
TwinTip-Holder for TransferMan NK for TransferMan	5177 208.010 5177 209.017
Microscope adapter see page 150	

